



# FLEXIBLE RELIABLE LIFECYCLE™

## LifeCycle™ Prenatal Screening software

Ideal solution for all prenatal screening laboratories.



### FLEXIBLE SCREENING

- Multiple screening strategies are supported
- Software is able to combine demographics and biomarker data to z-score results to produce risk



### TAILORED TO POPULATION

- Optimize risk calculation to your population
- Ability to level or change MoM equations



### ACCESSIBLE

- Accessible via web browser or over remote connection outside laboratory




### CONFIGURABLE


- High level of configurability and connectivity options



## LifeCycle™ software can produce z-score\* or risk



## Tailor risk calculation to your population with LifeCycle™ software



## LifeCycle™ software versions

Feature	v. 6.0	v. 7.0	v. 7.0 rev 2	v. 7.0 rev 3
User-configurable	✓	✓	✓	✓
Browser-based	✓	✓	✓	✓
Remote users	✓	✓	✓	✓
User-configurable file based LIMS import / export	✓	✓	✓	✓
T21, T18, T13 & pre-eclampsia risk	✓	✓	✓	✓
What if feature	✓	✓	✓	✓
Optional HL7 import / export for LIMS and Ultrasound equipment	✓	✓	✓	✓
Z-score calculation for Vanadis® NIPT system		✓	✓	✓
Fetal sex determination for Vanadis® NIPT system*		✓	✓	✓
Retest screening strategy for Vanadis® NIPT system			✓	✓
Ability to print labels for blood tubes			✓	✓
Pre-eclampsia twin risk calculation			✓	✓
sFlt-1/PIGF ratio calculation				✓

\*Sex determination is offered as an option for nonmedical purposes, only in countries where permitted. Sex determination is not offered in India under PNDT regulation.

**PerkinElmer, Inc.**  
 940 Winter Street  
 Waltham, MA 02451 USA  
 P: (800) 762-4000 or  
 (+1) 203-925-4602  
[www.perkinelmer.com](http://www.perkinelmer.com)

Wallac Oy  
 PO Box 10  
 20101 Turku, Finland  
 Phone: + 358 2 2678 111  
 Fax: + 358 2 2678 357

ISO 13485  
 ISO 14001  
 OHSAS 18001

Products may not be licensed in accordance with the laws in all countries, such as the United States. Please check with your local representative for availability.

For a complete listing of our global offices, visit [www.perkinelmer.com/ContactUs](http://www.perkinelmer.com/ContactUs)

Copyright © 2019–2021, PerkinElmer, Inc. All rights reserved. PerkinElmer® is a registered trademark of PerkinElmer, Inc. All other trademarks are the property of their respective owners.

1599-9802-06, February 2021 PKI